- **HANNA TREPINE, PUNCH BLOCK & ARTIFICIAL CHAMBER**
USER MANUAL1-10
- **TRÉPAN DE HANNA, BLOC « PUNCH » & CHAMBRE ARTIFICIELLE**
MANUEL UTILISATEUR 11-20
- **TRÉPANO DE HANNA, BLOQUE «PUNCH» & CÁMARA ARTIFICIAL**
MANUAL DEL USUARIO21-30
- **TRAPANO DI HANNA, BLOCCO « PUNCH » & CAMERA ARTIFICIALE**
MANUALE D'USO.....31-40
- **HANNA-TREPAN, STANZBLOCK & KÜNSTLICHE VORDERKAMMER**
BENUTZERHANDBUCH41-50
- **HANNA-TREPAAN, « PUNCH » BLOK & KUNSTMATIGE KAMER**
GEBRUIKERSHANDLEIDING51-60
- **TRÉPANO DE HANNA, BLOCO PERFURADOR & CÁMARA ARTIFICIAL**
MANUAL DO UTILIZADOR.....61-70
- **HANNA-TREPAN, STANSBLOCK & ARTIFICIELL KAMMARE**
ANVÄNDARMANUAL71-80
- **HANNA TREPANI, DONÖR « PANÇI » & YAPAY KAMARA**
KULLANICI KILAVUZU81-90





















Réf.: 65022

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









MORIA S.A. – 15, rue GEORGES BESSE - 92160 ANTONY – France
#65022 - D - 05/2017

SUMMARY

I. GENERAL RECOMMENDATIONS	2
II. HANNA TREPINE	
1. Description of the equipment.....	2
2. Technical specifications	2
3. Preparation before the operation.....	3
4. The operation	4
5. Troubleshooting	4
III. PUNCH BLOCK	
1. Description of the equipment.....	5
2. Technical specifications	5
3. Preparation before the operation.....	5
4. Trephination.....	6
5. Troubleshooting	6
IV. ARTIFICIAL CHAMBER	
1. Description of the equipment.....	7
2. Technical specifications	7
3. Preparation before the operation.....	8
4. Trephination.....	8
5. Troubleshooting	8
V. CLEANING AND STERILIZATION	9
VI. STORAGE & STERILIZATION BOX.....	9
VII. WARRANTY.....	10

I. GENERAL RECOMMENDATIONS

Optimal use of the Hanna Trepine, the punch block, and the artificial chamber requires strict compliance with the cleaning recommendations. If, for any reason, the equipment does not function properly, have it immediately checked by Moria. **We strongly recommend that the equipment be checked by Moria each year.**

The use of these devices with any single-use materials and/or components of a brand other than Moria's will immediately render the warranty null and void.

The latest version of the user manual and supplementary information for these systems is available on the Moria website: <http://www.moria-surgical.com>.

II. HANNA TREPINE

1. DESCRIPTION OF THE EQUIPMENT

Description	MORIA reference
Hanna Trepine (main body and suction cone)	17166
Disassembly key	17167
Depth adjustment key	17168
Storage and sterilization box	22517172
Hanna Trepine blade (diameter yyy mm) incision diameters available: 7.00mm / 7.25mm / 7.50mm / 7.75mm / 8.00mm / 8.25mm / 8.50mm / 8.75mm / 9.00mm / 9.50mm / 10.00mm / 10.50mm	17171Dyyyyx2 e.g.: the reference for a blade with an incision diameter of 7.00 mm is 17171D700x2

The Hanna Trepine is a high precision product. Its functionality and longevity are very dependent upon compliance with the instructions below.

The trephine consists of two parts: the main body with mechanism and the suction cone. It is supplied with an incision depth adjustment key (B) as well as a disassembly/assembly key (A) for use with blades over 8.5mm in diameter.

The unit is furnished with a storage and sterilization box in which each element has a specifically designated storage slot (see page 9). The artificial chamber can be stored and sterilized in this same box. **This box is not designed for transport and should be used only for storage and sterilization.**

In order to optimize the performance of the trephine, we advise you to use the dedicated artificial chamber.

2. TECHNICAL SPECIFICATIONS

	Dimensions	Weight	°C/storage	RH/storage
# 17166	Diam = 37 mm H = 32	60 g	Between 15°C/59°F and 35°C/95°F	Between 30 and 75%
# 22517172	L = 285 mm w = 125 mm H = 82 mm	1,4 Kg	Between 15°C/59°F and 35°C/95°F	Between 30 and 75%

3. PREPARATION BEFORE THE OPERATION

Ensure that the Hanna Trephine and its various elements are at room temperature before use.

ACTIONS	PRECAUTIONS
<p><u>Adjustment of the incision depth:</u></p> <ol style="list-style-type: none"> 1. Remove the suction cone (C) by unscrewing it from the main body of the trephine (D). 2. Check that the blade support (E) is in counter-clockwise position. <ol style="list-style-type: none"> a) Unlock the adjustment knob (F) by pulling gently (G), turn 90° so it is aligned and can fit between the rails (H), then leave the plate at the base of the adjustment knob in the setting position (I). b) Turn counter-clockwise until feeling resistance on the adjustment knob. You are now in «retracted» position. c) Lock the adjustment knob: reverse procedure in a) ⇒ realign the base of the adjustment knob between the rails, lift up and turn 90° again in the opposite direction in order to lock in upper position. 3. Use the key (B) to turn the ratchet collar (J). Turn clockwise to increase the incision depth (K): values from 0 to 1.2mm (that is 1200 µm). Note: the display (L) shows the figures from 0 to 12. 	<p>OPERATION TO BE CARRIED OUT BEFORE INSERTING THE BLADE (to avoid any risk of cuts and/or damaging the cutting edge of the blade)</p> <p>Never turn beyond number «12» (risk of significantly damaging the system).</p> <p>Setting the ratchet collar between two numbers, increases the incision depth by 50µm. Example: if the setting is between 5 and 6 ⇒ depth of 0.55mm, or 550 µm. The ratchet collar (J) cannot turn 360°. To return to «0», it must be turned counter-clockwise until the number «0» appears.</p>
<p><u>Inserting the blade :</u></p> <ol style="list-style-type: none"> 1. Turn the trephine blade locking ring counter-clockwise (M). 2. Insert the 3 small prongs on the blade into the 3 notches (N: blade unlocked). 3. Lock the blade by turning the ring clockwise (O: blade locked). 	<p>For blades with a diameter greater than 8.5mm, use the disassembly key (A) to remove the internal support ring (P) from the suction cone. (See photos of unscrewing internal support ring (Q) and the procedure. (R)). This enables visualization of the lowering of the larger blades.</p>
<p><u>Verification that the trephine is working correctly:</u></p> <ol style="list-style-type: none"> 1. Screw the suction cone (C) onto the main body of the trephine (D). 2. Attach the end of the silicon suction tubing to the hole in the suction cone (S). Fix the other end to a 50cc syringe (T). 3. Inspect the blade under a microscope. 4. Unlock the adjustment knob and turn it clockwise. The blade should move down to the chosen depth. When this depth is reached, you will feel resistance on the adjustment knob. It will continue to turn in order to complete the cut, but the trephine blade will not move down any further. 5. Withdraw the trephine blade completely (counter-clockwise) and lock the adjustment knob. The trephine is now ready for use. 	<p>Please ensure that the silicon tubing (U) is in good condition. We recommend replacing the tubing after each operation (references: #205 for one meter of silicon tube or #205B for a 15m reel of silicon tube).</p> <p>Any blade presenting a defect that may jeopardize the cutting performance must not be used.</p>

4. THE OPERATION

ACTIONS	PRECAUTIONS
<p>Positioning of the trephine:</p> <ol style="list-style-type: none"> 1. Check that the suction cone is firmly screwed onto the main body of the trephine. 2. Mark the center of the trephination either with an indentation or with a methylene blue surgical pen. Dry the cornea and the conjunctival cul-de-sac. The trephine is held in one hand while the other grasps the external rectus with tweezers to center the trephine guide on the corneal mark. 3. Press the trephine on the eye while the assistant applies suction with the syringe. Verify the vacuum of the trephine on the eye by pulling gently upwards. 	
<p>Trephination:</p> <ol style="list-style-type: none"> 1. Turn the adjustment knob clockwise. 2. When the Descemet's membrane is cut, a small discharge of aqueous humor appears (incomplete penetrations are possible by setting lower incision depths). Stop turning the adjustment knob and have the assistant release the suction to remove the trephine. 	<p>Other than the weight of the trephine itself, the surgeon must not exert any pressure on the eye during the operation.</p>
<p>End of the operation:</p> <ol style="list-style-type: none"> 1. After trephination is complete, withdraw the blade, turn the knob counter-clockwise, and lock it. 2. Unscrew the cone. Unlock the locking collar and remove the blade. Dispose of the blade in an appropriate container. 3. With the key (B), turn the ratchet collar (J) counter-clockwise in order to return to «0». 	<p>Warning! Failure to withdraw the blade will seriously damage the trephine and its blade depth regulation mechanism.</p> <p> Returning to the «0» position is compulsory (V). If not carried out, there is a high risk of jamming the system.</p>

5. TROUBLESHOOTING

- o If the device does not work properly, make sure there is no visible defect.
- o In the event of suction problems, check the integrity of the tubing and change it if necessary. Also check that the holes in the suction cone are not obstructed.
- o If the blade cannot be inserted properly, make sure that the guide has been removed (see precautions for inserting the blade page 3). If the blade depth cannot be set check for the presence of the guide.
- o If the device is locked, make sure that it has been returned to «0» using the ratchet collar (see. «end of the operation - precautions» in the table above).
- o In the case of cutting problems, check that the suction cone is securely screwed onto the main body of the trephine.
- o If the trephine is blocked, remove the blade completely from the trephine and lock the adjustment knob (turn the adjustment knob counter-clockwise until it stops, then turn it 90° clockwise and lock it). This process should be systematically carried out with the ratchet collar before each change in position.

III. PUNCH BLOCK

1. PUNCH BLOCK DESCRIPTION #17169

- ⇒ Teflon® stand (AA) with a Teflon® graft well (AB),
- ⇒ Cylindrical guide (AC),
- ⇒ 0.9mm (or 19 gauge) silicon tube (AD) + Luer connection,
- ⇒ Cutting system consisting of a cylinder (AE) and a locking ring (AF) at the end to secure the trephine blade.

2. TECHNICAL SPECIFICATIONS

	Dimensions	Weight	°C/storage	RH/storage
#17169 (mounted device)	L = 100 mm w = 100 mm H = 72 mm	514 g	Between 15°C/59°F and 35°C/95°F	Between 30 and 75%

3. PREPARATION OF THE PUNCH BEFORE THE OPERATION

Ensure that the punch base unit and its various elements are at room temperature before use. Ensure that the revolving nut is tight. The Teflon® well (AO) should be secure and even with the metal surface prior to punching.

ACTIONS	PRECAUTIONS
<ol style="list-style-type: none"> 1. Center the donor cornea, epithelial side down, on the Teflon® well. The epithelium is thus in contact with the Teflon® well. 2. Attach the silicon tube (Luer on the end - AG) to a 50cc syringe. Pull back on the plunger of the syringe to create a vacuum, flattening the donor cornea onto the well. 3. Inspect the single-use trephine blade under a microscope. 4. Mount the blade on the end (AF) : <ol style="list-style-type: none"> 4 -1 turn the locking ring to the open position (counter-clockwise) (AH), 4 -2 place the prongs of the blade onto the 3 markers (AK) (photo blade unlocked AI), 4 -3 lock the ring (clockwise). The blade is locked in position (AJ). 5. Screw the cylindrical guide (AC) onto the base (AA). 	<p>It is essential to check the condition of the surface of the well. It should be replaced if excessively scratched and at least after 10 uses.</p> <p>Ensure that the silicon tubing is in good condition. We recommend replacing it after each operation (references: #205 for one meter of silicon tube or #205B for a 15m reel of silicon tube).</p> <p>Any blade presenting a defect that may jeopardize the cutting performance must not be used.</p>

4. TREPHINATION

ACTIONS	PRECAUTIONS
1. Pull back on the plunger of the syringe to create a vacuum, causing the cornea to adhere to the Teflon® well.	If the adherence remains weak, it is possible to insert the trephine blade support (see step 2) and to apply suction again with the syringe. The weight of the trephine blade support enables adequate adherence.
2. Press the trephine blade support (AL) into the cylindrical guide (AC) so that the trephine blade rests on the endothelial surface of the cornea.	
3. Push the trephine blade support (AL) downward firmly with the fingertips in order to cut the graft (AM).	Do not use the palm of the hand; the resulting pressure may be excessive and could damage the Teflon® well. Teflon® enables the single-use blade to penetrate the material very slightly. This ensures that the graft is cut along the entire circumference.
4. Maintain the suction so that the graft remains firmly attached to the Teflon® well as you remove the trephine blade support (AL), then unscrew the cylindrical guide (AC).	

5. TROUBLESHOOTING

- o If the device does not function properly, make sure there is no visible defect. Check that the blade locking system slides correctly, by turning it in both directions.
- o In the case of a suction problem, check the integrity of the 2 seals and the tubing. Replace them if necessary.
- o In the case of an imprecise cut, check that the cylindrical guide is correctly screwed onto the base.
- o In case of particle deposits on the graft, check for the presence of Teflon® smears or shavings in the holes. Also check the condition of the Teflon® well and replace it if necessary.

Maintenance : Replacing the Teflon® well (ref. 17171P)

- Unscrew the revolving nut located under the base (AN).
- Push the well downward until it falls (AO).
- Place the new Teflon® well (as centrally as possible) without pressing on the nut (AP).
- Tighten the revolving nut. The Teflon® well slides into place automatically (AQ).

The Teflon® well must be replaced when damaged or after every 10 uses.

IV. ARTIFICIAL CHAMBER

WARNING: This chamber #18095C is specific to the Hanna Trephine. **It is not interchangeable with the «ALTK» #19161 artificial chamber.** Using the «ALTK» #19161 artificial chamber with the Hanna trephine can seriously damage the graft.

1. DESCRIPTION OF THE EQUIPMENT

The artificial chamber ref. #18095C can be used as an alternative to the Hanna punch block. When used with the Hanna trephine (ref. #17166), the chamber enables the trephination of the donor and recipient corneas on the epithelial side, while creating a consistent and symmetrical cut. This reduces disparities in shape and post-operative astigmatism. With the artificial chamber, it is possible to use the same size blade for the donor and the recipient, to cut the same side (epithelium), and to create a perpendicular cut which facilitates placing the graft on the patient.

The specific design of the trephine chamber stabilizes the cornea during trephination.

Note: The silicon irrigation tube (BE) is used in this technique not as a vacuum system, but as a system «under pressure» (approximately 50mmHg).

Please ensure that the silicon tubing is in good condition. We recommend replacing it after each operation (references: #205 for one meter of silicon tube or #205B for a 15m reel of silicon tube).

2. TECHNICAL SPECIFICATIONS

	Dimensions	Weight	°C/storage	RH/storage
#18095C	Diam = 73 mm H = 40 mm	326 g	Between 15°C/59°F and 35°C/95°F	Between 30 and 75%

3. PREPARATION OF THE ARTIFICIAL CHAMBER BEFORE THE OPERATION

ACTIONS	PRECAUTIONS
1. Insert a single-use blade of chosen diameter in the trephine (in certain cases it may be necessary to use an oversized blade - 0.25mm extra - for cutting the graft).	
2. Select the maximum incision depth (1.2mm) on the trephine for complete perforation.	
3. Remove the retaining ring (BA) by loosening the locking collar (BB) (clockwise) in the artificial chamber (BF =loosening of the ring; BG = chamber without ring).	
4. Apply visco-elastic to the graft endothelium. The graft is then centered on the support (BC), endothelial side facing downwards.	Before placing the graft on its support, it is recommended to fill the chamber with liquid (e.g.: BSS™ - registered trademark from Alcon® or saline solution) in order to leave as little air as possible during step no. 7.
5. Cover the graft support (BC) by installing the retaining ring (BA) on the support block (BD).	
6. Lock the retaining ring on the graft (BA) using the locking collar (BB) – counter-clockwise. The ring (BA) must be held when closing the locking collar (BB).	
7. Fill the chamber with liquid (such as BSS™ or saline solution) through the irrigation tube in order to obtain adequate pressure (approximately 50mmHg). Expel the air bubbles under the graft by opening and then closing the purge tap (BH) on the retaining ring (BA).	

4. TREPHINATION

ACTIONS	PRECAUTIONS
1. Place the pre-calibrated Hanna trephine (BI) on the retaining ring (BA) in the artificial chamber (BJ), provided for this purpose.	Check that the trephine is correctly pressed into the retaining ring (metal against metal). There is no need to apply suction during incision.
2. Cut the graft (in the same way as the recipient' cornea.	
3. Remove the trephine from its support.	
4. Stop injecting the liquid.	
5. Grasp the retaining ring (BA) and unlock the ring (BB).	
6. Remove the retaining ring (BA) from its support.	
7. Remove the graft and protect it (before surgery).	

5. TROUBLESHOOTING

In the event of suction problems:

- o Check the integrity of the seal, as well as the correct functioning of the purge screw and its seal.
- o Check that the retaining ring is correctly maintained by the locking collar.

V. CLEANING AND STERILIZATION

The recommendations for cleaning, disinfection and sterilization are updated and available on Moria's website: <http://www.moria-surgical.com>.

1. Disassembly: The device must be completely disassembled before decontamination and cleaning operations.
2. Disinfection: Immerse the individual parts as quickly as possible after surgery in a decontaminating solution (pre-disinfectant alkaline detergent) for at least 15 min.
3. Rinsing: Rinse the device thoroughly with water below 30°C/86°F.
4. Cleaning: Use an ultrasound tank. Place the various parts in separate plastic tanks so as not to bring them in direct contact with the walls of the tank. Leave them in the tanks for approximately 10 minutes.
5. Drying: Whichever cleaning method is used, the instruments must be perfectly dry before further use, in order to avoid any recontamination.
6. Sterilization: The device can be autoclave sterilized cycle up to 134°C (273°F)/18 min.
7. Storage: Store in an appropriate, dry place.

For optimal cleaning of the trephine: Clean the trephine and the suction hole with distilled water at a temperature below 30°C/86°F. Ensure that no debris remains. Disassemble the black part of the guide (W) with the disassembly key (A) - see photo of unscrewing (X). Then clean the various parts with a soft toothbrush.

For optimal cleaning of the punch: Unscrew the nut located under the base in order to separate the Teflon® disc and the silicon tube before any cleaning operations.

 **Note concerning NCTA:** for France, follow the recommendations of the Ministry of Health concerning the deactivation of NCTA (DGS/5 C/DHOS/E 2 n° 2001-138 du 14/03/01). Moria does not recommend chemical inactivation using BLEACH. Only inactivation using CAUSTIC SODA is recommended for rustproof materials and titanium. The punch block as well as the artificial chamber withstand inactivation using caustic soda. For the Hanna trephine, only the suction cone withstands inactivation using caustic soda.

In the case of decontamination using caustic soda, **ensure thorough and prolonged rinsing of the device with distilled water. Insufficient rinsing (remains of caustic soda) can seriously damage the cornea.**

VI. STORAGE AND STERILIZATION BOX

The storage and sterilization box for the Hanna Trephine, the punch block, and the artificial chamber is the reference **22517172**:

- photo BK of the bottom of the empty box #22517172
- photo BL of the internal elements - punch version
- photo BM of the internal elements - punch version - Teflon® cover on the trephine (**warning: do not forget Teflon® protection on the base of the trephine!**)
- photo BN of the internal elements - artificial chamber version.

VII. WARRANTY

Scope of application of the warranty

Designation	MORIA references
Hanna Trephine (main body and suction cone)	17166
Disassembly key	17167
Incision depth adjustment key	17168
Punch block	17169
Artificial chamber	18095C
Storage and sterilization box	22517172

- The above items as well as spare parts and labor necessary for their repair are covered by warranty. Any items returned must be sent in their original packaging, **after having previously been decontaminated**.
- The maintenance operations and the replacement of spare parts will be exclusively carried out by technicians authorized by MORIA.
- Instruments and accessories that cannot be reused such as the blades and the Teflon® well, are excluded from this warranty.

Case of non-application of the warranty

- Defects or malfunction that occur out of the warranty period (defined below).
- Normal wear and tear.
- Negligence or usage that does not comply with the specifications in the user manual.
- The use of supplies, spare parts, or accessories other than those supplied by MORIA (for example: blades and tubing not supplied by MORIA).
- Any disassembly, modification, or intervention carried out on the devices by a person not authorized by MORIA.

Warranty period

- The warranty takes effect on the date the material is dispatched.
- The duration of the warranty is 12 months from the date of effect.

Liability

- The liability of MORIA is limited to the supply of the services mentioned in paragraph «scope of application of the warranty». MORIA will not be held responsible for any direct or indirect damage notably financial, suffered by the client owing to the interventions within the scope of this warranty.
- For any dispute concerning the interpretation or the execution of the present contract or the present general terms and conditions, the Commercial Court of Nanterre, France, will have sole jurisdiction.

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